

# Colostomy Bag

C1NB10-30 and C1NB10V-HC-30




## Skin Wafer

The skin friendly and very flexible skin wafer is made of 100% hydrocolloid. It provides a secure fit, protects the skin & absorbs moisture efficiently. There is excellent balance between secure adhesion and comfortable removal.



## Tissue

The very soft, breathable and water resistant tissue is comfortable against the skin and prevents perspiration forming underneath the pouch. It dries very quickly after swimming or showering.



## Filter

Our exceptionally good activated carbon filter provides a very long odour neutralising pathway and can cleanse large amounts of gas. This means the gas will be odour free when it leaves the filter.



## Viewing Window

The ostomy pouch includes a viewing window which makes it possible to inspect the stoma as well as the output from the stoma. It is also very helpful when applying the pouch to make sure it is located correctly over the stoma.



## Integrated Vent

The integrated vent is a new and novel way to vent the pouch allowing air out or in, for instance in the case of a vacuum occurring. It can also be used for more effective cleansing by rinsing the pouch via the vent opening.

## Adhesive Strength

- Comfortable skin wafer with excellent adhesion for secure and guaranteed adhesion
- Extra secure adhesion, even on the edges

## Filter and Air Vent

- Unique integrated vent to control the air
- Durable filter control of odours

## Skin Comfort

- The comparative of skin-friendliness: skin comfort
- Extremely flexible and comfortable