



心电电缆及附件 ECG Cables & Accessories 2014

Content

1. ECG Monitoring Cables & Lead wires

1.1. ECG Cables & Lead wires Cross reference

1.2. Major types or brands ECG lead wires

2. Telemetry ECG Cables

3. Disposable ECG lead wires

4. ECG Lead Adapter

5. Diagnosis ECG Cables & Accessories

5.1. One piece EKG fixed cable

5.2. EKG Trunk Cable

5.3 EKG Lead wire

5.4 EKG Accessories

Compliance with 符合:

ANSI/AAMI EC53-1998(R)2001

EN 60601-1:1990

ISO10993-5:2003E



FDA 510(K) clearance No:K082959

Patient Cable & Leadwire 病人电缆&导联线



SFDA 医疗器械注册证

证书号码 No.: 1210080

<https://endo.id>

内容

1. 心电监护电缆与导联线

1.1 心电电缆导联线对照表

1.2 主要类型或品牌心电导联线

2. 心电遥测电缆

3. 一次性心电导联线

4. 心电导联转换头

5. 诊断用心电电缆及附件

5.1 一体化心电电缆

5.2 心电主电缆

5.3 心电导联线

5.4 附件

1.ECG Monitoring Cables & Lead wires 心电监护电缆与导联线

Compatibility 兼容性:

Major brands ECG or multi-parameter monitoring system

主要品牌的心电和多参数监护系统
















Specification 规格

Cable material 线材	TPU Jacket, Light gray TPU 外壳, 浅灰色	Connector &Yoke 接头和分线器	Nylon plug & Golden plated pins 尼龙插头和镀金插针
Main Cable length 主电缆长度	2.5m	Lead wire length 导联线长度	0.9m

Shielded Cable & lead wire 屏蔽型电缆和导联线	Figures 图片	REF# examples 型号	Reference to REF# 型号命名参考规则
One piece Fixed Cable 一体化电缆 (Shielded 屏蔽)		2586P	2: Shielded(1=unshielded) 2: 屏蔽 (1=非屏蔽) 5: 5-lead(3= 3 lead) 5: 5-导联(3=3-导联) 86: factory/model code 86: 原厂/规格型号 N: with neonate clip end(Default is adult) N: 新生儿夹(缺省: 成人) See the next page for the figure(见下页中的图片) S: snap end(P=grabber end) S: 扣式末端(P=钳式) I: IEC code (Default is AHA) I: 欧标 (缺省: 美标)
Trunk Cable 主电缆 (Shielded 屏蔽)		AA-2585	AA: cable yoke type, up to lead wire AA: 分线盒类型,取决于导联线 2: Shielded(1=unshielded) 2: 屏蔽 (1=非屏蔽) 5:5=lead(3= 3 lead) 5: 5-导联(3=3-导联) 85: factory/model code 86: 原厂/规格型号 Default : AHA code 缺省: 美标
Lead wire 导联线 (Shielded 屏蔽)		NKB6-90S-I	NK: lead wire type(Acc to yoke) NK: 导联线类型 (按分线盒) 6: 6-lead (3= 3 lead,5=5lead) 5: 5-导联(3=3-导联) 90: 90cm(wire length) 90: 90cm(导联线长度) S: snap end(P=grabber end) S: 扣式末端 (P=钳式) I: IEC code(Default : AHA code) I: 欧标 (缺省: 美标)

Unshielded Cable & lead wire 非屏蔽型电缆和导联线	Photo 图片	REF# 型号	Reference to REF# 型号命名参考规则
Din trunk Cable D 型主电缆 (1.5mm, Unshielded 非屏蔽) Used with Din lead wire Especially for neonate ECG monitoring 与 D 型导联线配合使用，常用在新生儿 ECG 监护系统		D-1540-I	D: Din style cable 1: Unshielded 5: 5-lead(3= 3 lead) 40: factory/model code I: IEC code (Others see below) D: D 型电缆 1: 非屏蔽 5: 5-导联(3=3-导联) 40: 原厂/规格型号 I: 欧标 (其它如下)
Din lead wire D 型导联线 (1.5mm, Unshielded 非屏蔽) Used with Din trunk cable only 仅与 D 型主电缆配合使用		D5-90S-I	D: Din lead wire (Acc to yoke) 5: 5-lead (3= 3 lead) 90: 90cm(wire length) S: snap end(P=grabber end) I: IEC code (Others see below) D: D 型导联线 (按分线盒) 5: 5-导联(3=3-导联) 90: 90cm(导联线长度) S: 扣式末端 (P=钳式) I: 欧标 (其它如下)
Neonate lead wire 新生儿导联线 (1.5mm, Unshielded 非屏蔽) Used with Din trunk cable only 仅与 D 型主电缆配合使用		D3-90NP-I	D: Din lead wire (Acc to yoke) 3: 3-lead (5= 5 lead) 90: 90cm(wire length) NP=mini grabber end only I: IEC code (Others see below) D: D 型导联线 3: 3-导联(5=5-导联) 90: 90cm(导联线长度) NP=迷你钳式末端 I: 欧标 (其它如下)

ECG Color Code on electrodes (Regular grabber/Mini grabber/Regular snap)			ECG 电极端颜色规则 (常规钳式/迷你钳式/常规扣式)		
AHA(A)	Alternate AHA(H)	Standard IEC(I)	Alternate IEC(IH)	AHA(A)	Standard IEC(I)
					
					
					

1.1. ECG Cables & Lead wires Cross reference 心电电缆导联线对照表

OEM &Models OEM 原厂型号	Monitor Plug 仪器端接头	Lead No. 导联数	Trunk Cable 主电缆	Lead yoke Conn. 导联线接头	Lead wire/w Snap 导联线：扣式	Lead wire/w Grabber 导联线：钳式	One piece/w Snap 一体化扣式	One piece/w Grabber 一体化钳式	Din Trunk Cable D 型主电缆
AAMI									
Without resistance BCI Datascope MDE/Invivo Mindray Spacelabs WelchAllyn Protocol		3	D-1340		D3-90S	D3-90P	2340S	2340P	
			D-1340-I		D3-90S-I	D3-90P-I	2340S-I	2340P-I	
		Price	C1		A8	B3	C1	C1	
		5	D-1540		D5-90S	D5-90P	2540S	2540P	
			D-1540-I		D5-90S-I	D5-90P-I	2540S-I	2540P-I	
		Price	C2		B4	B6	C2	C2	
With 1K resistance Same as above		3	D-1340R		D3-90S	D3-90P	2340RS	2340RP	
			D-1340R-I		D3-90S-I	D3-90P-I	2340RS-I	2340RP-I	
		Price	C1		A8	B3	C1	C1	
		5	D-1540R		D5-90S	D5-90P	2540RS	2540RP	
			D-1540R-I		D5-90S-I	D5-90P-I	2540RS-I	2540RP-I	
		Price	C2		B4	B6	C2	C2	
LL style		3	LL-2340		LL3-90S	LL3-90P	2340S	2340P	D-1340
			LL-2340-I		LL3-90S-I	LL3-90P-I	2340S-I	2340P-I	D-1340-I
		Price	C6		B7	B7	C1	C1	C1
		5	LL-2540		LL5-90S	LL5-90P	2540S	2540P	D-1540
LL-2540-I	LL5-90S-I		LL5-90P-I		2540S-I	2540P-I	D-1540-I		

PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

		Price	C6		C2	C2	C2	C2	C2
OEM &Models OEM 原厂型号	Monitor Plug 仪器端接头	Lead No. 导联数	Trunk Cable 主电缆	Lead yoke Conn. 导联线接头	Lead wire/w Snap 导联线: 扣式	Lead wire/w Grabber 导联线: 钳式	One piece/w Snap 一体化扣式	One piece/w Grabber 一体化钳式	Din Trunk Cable D 型主电缆
Artema/S&W									
8000,Diascope1,-2, Diascope NT, Athena		3	NA	NA	NA	NA	2397S	2397P	D-1397
		Price					2397S-IH	2397P-IH	D-1397-IH
							C4	C4	C4
BCI	See AAMI								
Biolight 宝莱特 (China)									
M7000,M9500,AS Modular		5	NA	NA	NA	NA	2557S	2557P	D-1557
		Price					2557S-I	2557P-I	D-2557-I
							C6	C6	C6
Biolight AnyView A2/A3/A5/A6/A8 Q3/Q5/Q7 V6 Vital Sign Monitor		3	NA	NA	NA	NA	2367S	2367P	D-1367
		Price					2367S-I	2367P-I	D-1367-I
Bionet(Korea)									
		5	NA	NA	NA	NA	2514S	2514P	D-1514
		Price					2514S-I	2514P-I	D-2514-I
							C6	C6	C6
Biosys(Korea)									
		3	NA	NA	NA	NA	2335S	2335P	D-1335
		price					2335S-I	2335P-I	D-1335-I
							C4	C4	C4

PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

CAS MEDICAL									
		3	NA	NA	NA	NA	2339S	2339P	D-1339
							2339S-IH	2339P-IH	D-1339-IH
		price					C4	C4	C4
OEM &Models OEM 原厂型号	Monitor Plug 仪器端接头	Lead No. 导联数	Trunk Cable 主电缆	Lead yoke Conn. 导联线接头	Lead wire/w Snap 导联线：扣式	Lead wire/w Grabber 导联线：钳式	One piece/w Snap 一体化扣式	One piece/w Grabber 一体化钳式	Din Trunk Cable D 型主电缆
Colin/Omron(Japan)									
BP 308/308E/408	See AAMI								
BP88S		3	CL-2520		CL3-90S	CL3-90P	2320S	2320P	D-1320
			CL-2520-I		CL3-90S-I	CL3-90P-I	2320S-I	2320P-I	D-1320-I
		Price	C6		B7	B7	C4	C4	C4
		5	CL-2520		CL5-90S	CL5-90P	2520S	2520P	D-1520
			CL-2520-I		CL5-90S-I	CL5-90P-I	2520S-I	2520P-I	D-1520-I
	Price	C6		C2	C2	C5	C5	C5	
BP88 BP306		3	CL-2521		CL3-90S	CL3-90P	2321S	2321P	D-1321
			CL-2521-I		CL3-90S-I	CL3-90P-I	2321S-I	2321P-I	D-1321-I
		Price	C6		B7	B7	C4	C4	C4
		5	CL-2521		CL5-90S	CL5-90P	2521S	2521P	D-1521
			CL-2521-I		CL5-90S-I	CL5-90P-I	2521S-I	2521P-I	D-1521-I
	Price	C6		C2	C2	C5	C5	C5	
Criticare	See AAMI								
602A-A, Poet 2, 504US, Scholar,506,507,508, 1100,1508,8100									
Creative 科瑞康(China)									
		3					2330S	2330P	D-1330
							2330S-I	2330P-I	D-1330-I


PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

		Price					C1	C1	C1
		5					2530S	2530P	D-1530
							2530S-I	2530P-I	D-1530-I
		Price					C2	C2	C2
OEM &Models OEM 原厂型号	Monitor Plug 仪器端接头	Lead No. 导联数	Trunk Cable 主电缆	Lead yoke Conn. 导联线接头	Lead wire/w Snap 导联线: 扣式	Lead wire/w Grabber 导联线: 钳式	One piece/w Snap 一体化扣式	One piece/w Grabber 一体化钳式	Din Trunk Cable D 型主电缆
Datascope /Mindray	See AAMI								
Datascope -Mindray V12, V21		3	DT-2526 0012-00-1745-01		DT3-90S	DT3-90P	2326S	2326P	D-1326
			DT-2526-I		DT3-90S-I	DT3-90P-I	2326S-I	2326P-I	D-1326-I
		Price	C6		C2	C2	C4	C4	C4
		5	DT-2526		DT5-90S	DT5-90P	2526S	2526P	D-1526
			DT-2526-I		DT5-90S-I	DT5-90P-I	2526S-I	2526P-I	D-1526-I
		Price	C6		C2	C2	C6	C6	C6
Dixtal									
EP12			NA	NA	NA	NA	21070S	21070P	NA
			Na	NA	NA	NA	21070S-I	21070P-I	NA
		Price					D3	D3	
Drager									
UM1/3/3.1, RM1, PM8010, P M8014, Babyguard8000, C ato, PB8800, Series 8000, Cicero, Parameterbox		3	NA	NA	NA	NA	2377S	2377P	D-1377
							2377S-I	2377P-I	D-1377-I
		Price					C4	C4	C4
ECAFIX									
		5	NA		NA	NA	2552S 2552S-I	2552P 2552P-I	NA

PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

							C2	C2	
Fukuda Danshi									
		5	SM-2599		SM5-90S	SM5-90S	2599S	2599P	D-1599
			SM-2599-I		SM5-90S-I	SM5-90S-I	2599S-I	2599P-I	D-1599-I
		Price	C6		C2	C2	C6	C6	C6
		12	MQ-21099		MQ5-90S	MQ5-90P	21099S	21099P	
			MQ-21099-I		MV5-90S	MV5-90P			
			MQ5-90S-I		MQ5-90P-I	21099S-I	21099P-I		
		Price	D3		2XC2	2XC2	D3	D3	
DS5100E , DS5300W , DS5400、DS7001 Mindray/Datascope /Expert		5					2559S	2559P	D-1559
							2559S-I	2559P-I	D-1559-I
		Price				C6	C6	C6	
GE/Critikon									
Dianmap Pro 1000 MPS Select Portable monitors		3	DX-2340R		DX3-90S	DX3-90P	2340RS	2340RP	D-1340R
			DX-2340R-I		DX3-90S-I	DX3-90P-I	2340RS-I	2340RP-I	D-1340R-I
		Price	C6		B7	B7	C1	C1	C1
GE-Datex-Ohmeda									
AS/3, CS/3, Cardiocap II, Satlite Plus S5 Without Resister for respiration		3	DX-2392		DX3-90S	DX3-90P	2392S	2392P	D-1392
			DX-2392-I		DX3-90S-I	DX3-90P-I	2392S-I	2392P-I	D-1392-I
		Price	C6		B7	B7	C4	C4	C4
		5	DX-2592		DX5-90S	DX5-90P	2592S	2592P	D-1592
			DX-2592-I		DX5-90S-I	DX5-90P-I	2592S-I	2592P-I	D-2592-I
		Price	C6		C2	C2	C6	C6	C6
AS/3, CS/3, Cardiocap II, Satlite Plus S5		3	DX-2395		DX3-90S	DX3-90P	2395S	2395P	D-1395
			DX-2395-I		DX3-90S-I	DX3-90P-I	2395S-I	2395P-I	D-1395-I
		Price	C6		B7	B7	C4	C4	C4

PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

With resisters		5	DX-2595		DX5-90S	DX5-90P	2595S	2595P	D-1595
			DX-2595-I		DX5-90S-I	DX5-90P-I	2595S-I	2595P-I	D-2595-I
		Price	C6		C2	C2	C6	C6	C6
		12	DX-21095		MQ5-90S	MQ5-90P	21095S	21095P	
			DX-21095-I		MQ5-90S-I	MQ5-90P-I	21095S-I	21095P-I	
		Price	NA		C2X2	C2X2	NA	NA	
GE/Hellige									
Eagle 1000 Def.. SCP 900/ Cardioserv		3	DX-2355		DX3-90S	DX3-90P	2355S	2355P	D-1355
			DX-2355-I		DX3-90S-I	DX3-90P-I	2355S-I	2355P-I	D-1355-I
		Price	C6		B7	B7	C4	C4	C4
		5	DX-2555		DX5-90S	DX5-90P	2555S	2555P	D-1555
			DX-2555-I		DX5-90S-I	DX5-90P-I	2555S-I	2555P-I	D-1555-I
Price	C6	C2	C2	C6	C6	C6			
SMC,SMS,SMK,SMU,SMV									
TRAM 7200		5	MQ-2584		MQ5-90S	MQ5-90P	2584S	2584P	D-1584
			MQ-2584-I		MQ5-90S-I	MQ5-90P-I	2584S-I	2584P-I	D-1584-I
		price	C6		C2	C2	C6	C6	C6
GE/Marquette									
Eagle,Solar,Dash monitor MAC-Lab Cath Lab Datex-Ohmeda M-Prestn & “E” modules	 	3	MQ-2586		MQ3-90S	MQ3-90P	2386S	2386P	D-1386
			MQ-2586-I		MQ3-90S-I	MQ3-90P-I	2386S-I	2386P-I	D-1386-I
		Price	C6		C2	C2	C4	C4	C4
		5	MQ-2586		MQ5-90S	MQ5-90P	2586S	2586P	D-1586
			MQ-2586-I		MQ5-90S-I	MQ5-90P-I	2586S-I	2586P-I	D-1586-I
		Price	C6		C2	C2	C6	C6	C6
		6	MQ-2686		MQB6-90S	MQB6-90P	NA	NA	NA
	MQ-2686-I	MQB6-90S-I	MQB6-90P-I		NA	NA	NA		

<https://end>


PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

		Price	C8.2		C4	C4	D3	D3	
		12	MQ-21086		MQ5-90S MV5-90S	MQ5-90P MV5-90P	MQ-21086S	MQ-21086P	
			MQ-21086-I		MQ5-90S-I MV5-90S-I	MQ5-90P-I MV5-90P-I	MQ-21086S-I	MQ-21086P-I	
		Price	D3		C2X2	C2X2	D3	D3	
Kontron									
Micromon, Minimon Vitalmon, Supermon Kolormon		3	NA		NA	NA	2398S	2398P	D-1398
							2398S-IH	2398P-IH	D-1398-IH
		Price					C4	C4	C4
105, 107, 108, 125, 128, VSM510。 Old generation		3	NA		NA	NA	2393S	2393P	D-1393
							2393S-I	2393P-IH	D-1393-IH
		Price					C4	C4	C4
Lohmerier									
M010,M011/M211 Defi 501		3			LM3-90S	LM3-90P	2358S	2358P	D-1358
					LM3-90S-I	LM3-90P-I	2358S-I	2358P-I	D-1358-I
		Price			B7	B7	C4	C4	C4
		5			LM5-90S	LM5-90P			
					LM5-90S-I	LM5-90P-I			
		Price			C2	C2			
MEK									
MP500, MP600, MP1000		3	NA		NA	NA	2310S	2310P	D-1310
							2310S-I	2310P-I	D-1310-I
		Price					C2	C2	C2
		5	NA		NA	NA	2510S	2510P	D-1510
								2510S-I	2510P-I


PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia



		Price					C4	C4	C4
		3	NA		NA	NA	2311S	2311P	D-1311
							2311S-I	2311P-I	D-1311-I
		Price					C1	C1	C1
		5	NA		NA	NA	2511S	2511P	D-1511
							2511S-I	2511P-I	D-1511-I
		Price					C2	C2	C2
Mennen									
Series 9000,1100,2000		3	N/A		NA	NA	2374S	2374P	D-1374
							2374S-I	2374P-I	D-1374-I
		Price					C4	C4	C4
		5	N/A		NA	NA	2574S	2574P	D-1574
							2574S-I	2574P-I	D-1574-I
		Price					C6	C6	C6
Envoy		3	AA-2375		MN3-90S	MN3-90P	2375S	2375P	D-1375
			AA-2375-I		MN3-90S-I	MN3-90P-I	2375S-I	2375P-I	D-1375-I
		Price	C6		B7	B7	C4	C4	C4
		5	AA-2575		MN5-90S	MN5-90P	2575S	2575P	D-1575
			AA-2575-I		MN5-90S	MN5-90P-I	2575S-I	2575P-I	D-1575-I
		Price	C6		C2	C2	C6	C6	C6
Mindray									
PM7000/8000/9000	See AAMI								
PM5000		3	MR-2512		MR3-90S	MR3-90P	2312S	2312P	D-1312
PM6000			MR-2512-I		MR3-90S-I	MR3-90P-I	2312S-I	2312P-I	D-1312-I
BeneViewT5/T6/T8		Price	C6		B7	B7	C4	C4	C4

PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

T5/T6/T8		5	MR-2512		MR5-90S	MR5-90P	2512S	2512P	D-1512
			MR-2512-I		MR5-90S-I	MR5-90P-I	2512S-I	2512P-I	D-1512-I
		Price	C6		C2	C2	C6	C6	C6
NEC									
 		3			NECA3-90S	NECA3-90P			
					NECA3-90S-I	NECA3-90P-I			
		Price		B7	B7				
		3			NECB3-90S	NECB3-90P			
					NECB3-90S-I	NECB3-90P-I			
		Price		B7	B7				
		3			NECC3-90S	NECC3-90P			
					NECC3-90S-I	NECC3-90P-I			
		Price			B7	B7			
Nihon Kohden									
Def. Cardiolife Lifescope series 6-8-9-11 6102-6205-6120		3	NK-2508		NKA3-90S	NKA3-90P	2308S	2308P	D-1308
			NK-2508-I		NKA3-90S-I	NKA3-90P-I	2308S-I	2308P-I	D-1308-I
		Price	C6		B7	B7	C4	C4	C4
		5	NK-2508		NKA5-90S	NKA5-90P	2508S	2508P	D-1508
			NK-2508-I		NKA5-90S-I	NKA5-90P-I	2508S-I	2508P-I	D-1508-I
		Price	C6		C2	C2	C6	C6	C6
		3	NK-2507		NKA3-90S	NKA3-90P	2307S	2307P	D-1307
			NK-2507-I		NKA3-90S-I	NKA3-90P-I	2307S-I	2307P-I	D-1307-I
		Price	C6		B7	B7	C4	C4	C4
		5	NK-2507		NKA5-90S	NKA5-90P	2507S	2507P	D-1507
			NK-2507-I		NKA5-90S-I	NKA5-90P-I	2507S-I	2507P-I	D-1507-I
		Price	C6		C2	C2	C6	C6	C6

PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia


TEC-7621C/TEC-7631C TEC-7721C/TEC-7731C 7700		3					2313S	2313P	
							2313S-I	2313P-I	
		Price					C6	C6	
BSM-2300 lifescope I BSM-4100 lifescope P BSM-5100 lifescope A BSM-9510 lifescope M BSM-9800 lifescope S		3	NK-2509		NKB3-90S	NKB3-90P	2309S	2309P	D-1309
			NK-2509-I		NKB3-90S-I	NKB3-90P-I	2309S-1	2309P-I	D-1309-I
		Price	C8		B7	B7	C6	C6	C6
		6	NK-2509		NKB6-90S	NKB6-90P	2609S	2609P	
			NK-2509-I		NKB6-90S-I	NKB6-90P-I	2609S-1	2609P-I	
		Price	C8		C3	C3	C8	C8	
Life Scope: OPV-1500		3	(JC-103TA)		NKB3-90S	NKB3-90P	2316S	2316P	D-1316
			(JC-103T)		NKB3-90S-I	NKB3-90P-I	2316S-1	2316P-I	D-1316-I
		Price	NA		B7	B7	C6	C6	C6
		3	NK-2540		NKA3-90S	NKA3-90P	2340S	2340P	D-1340
			NK-2540-I		NKA3-90S-I	NKA3-90P-I	2340S-I	2340P-I	D-1340-I
		Price	C6		B7	B7	C1	C1	C1
		5	NK-2540		NKA5-90S	NKA5-90P	2540S	2540P	D-1540
			NK-2540-I		NKA5-90S-I	NKA5-90P-I	2540S-I	2540P-I	D-1540-I
		Price	C6		C2	C2	C2	C2	C2
		Philips/HP							
43100,43200A, 77020A/AC/AR		3	AA-23029		AA3-90S	AA3-90P	23029S	23029P	D-13029
			AA-23029-I		AA3-90S-I	AA3-90P-I	23029S-I	23029P-I	D-13029-I

Serie 78300, serie 78500 Serie 78600, M2400A,M2406A, M2424A		Price	C6		B7	B7	C4	C4	C4
		5	AA-25029		AA5-90S	AA5-90P	25029S	25029P	D-15029
			AA-25029-I		AA5-90S-I	AA5-90P-I	25029S-I	25029P-I	D-15029-I
		Price	C6		C2	C2	C6	C6	C6
Merlin 78352C/78354C/ 78834CM1001A/ M1002A, Virdia, Ominicare CMA, ML65A,66A Heartstream Codemaster XL ...		3	AA-2385		AA3-90S	AA3-90P	2385S	2385P	D-1385
			AA-2385-I		AA3-90S-I	AA3-90P-I	2385S-I	2385P-I	D-1385-I
		Price	C6		B7	B7	C4	C4	C4
		5	AA-2585		AA5-90S	AA5-90P	2585S	2585P	D-1585
			AA-2585-I		AA5-90S-I	AA5-90P-I	2585S-I	2585P-I	D-1585-I
		Price	C6		C2	C2	C6	C6	C6
Modules: M1001B, M1002B, Multi- Measurement Modules: M3000A, M3001A, M3002A, MP2: M8102A, MP5: M8105A, SureSigns Monitors: 863064, 863065, 863066, 863068, 863087, 863088		3	HP-2385		HP3-90S	HP3-90P	2385S	2385P	D-1385
			HP-2385-I		HP3-90S-I	HP3-90P-I	2385S-I	2385P-I	D-1385-I
		Price	C7		C2	C2	C4	C4	C4
		5	HP-2585		HPA5-90S	HPA5-90P	2585S	2585P	D-1585
			HP-2585-I		HPA5-90S-I	HPA5-90P-I	2585S-I	2585P-I	D-1585-I
		Price	C8		C4	C4	C6	C6	C6
Multi-Measurement Modules: M3000A,		10	HP-21085		HPA5-90S(blue) HPB5-90S(Red)	HPA5-90P(Blue) HPB5-90P(Red)	HP-21085S	HP-21085P	



er, distributor, and manufacturer of medical products in Indonesia

C18; Patient Monitors: MP5: M8105A Opt C12, MP2: M8102A Opt C12, ECG Modules: M1001B, M1002B, 860335, 860338			HP-21085-I		HPA5-90S-I(Blue) HPB5-90S-I(Red)	HPA5-90P-I(Blue) HPB5-90P-I(Red)	HP-21085S-I	HP-21085P-I	
	Price	D3			C4X2	C4X2	D3	D3	
Physio Control									
Lifepack5,6,7,8,9,9P,VSM 2,1440,Vingmed 5	See GE Critikon								
LIFEPAK12、LIFEPAK20E、 LIFEPAK15 With resisters 		3					2327S	2327P	
							2327S-I	2327P-I	
		Price					C4	C4	
		5					2527S	2527P	
							2527S-I	2527P-I	
		Price					C6	C6	
		12	MD-21027		MDA4-120S MDB6-90S		21027S	21027P	
			MD-21027-I		MDA4-120S-I MDB6-90S-I		21027S-I	21027P-I	

		Price	D3		D3		D3	D3	
Schiller	NA								
Siemens									
Sirecust series 400 Sirecust 341,610/620/630/ 720/730/732/900/960/961 1260/1281/1801, 400-500 ECG		3	SM-2376		SM3-90S	SM3-90P	2376S	2376P	D-1376
			SM-2376-I		SM3-90S-I	SM3-90P-I	2376S-I	2376P-I	D-1376-I
		Price	C6		B7	B7	C4	C4	C4
		5	SM-2576		SM5-90S	SM5-90P	2576S	2576P	D-1576
			SM-2576-I		SM5-90S-I	SM5-90P-I	2576S-I	2576P-I	D-1576-I
		Price	C6		C2	C2	C6	C6	C6
SC5000,SC6000, SC7000 SC8000,SC9000 Multi-Parameter Cable SM-STE5 		3	SM-STE5	See the following page	SMB3-90S	SMB3-90P			
			SM-STE5		SMB3-90S-I	SMB3-90P-I			
		Price	D8		B7	B7			
		5	SM-STE5		SM5-90S	SM5-90P			
			SM-STE5		SM5-90S-I	SM5-90P-I			
		Price	D8		C2	C2			

PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

		3	SM-2323		SM3-90S	SM3-90P	2323S	2323P	D-1323
			SM-2323-I		SM3-90S-I	SM3-90P-I	2323S-I	2323P-I	D-1323-I
		Price	C8		B7	B7	C6	C6	C6
Spacelabs									
90701,90303,90305,90308,,90309,90385	See AAMI 								
Ultraview 90492, 91220, 90496, 91496		3	SL-2396		SL3-90S	SL3-90P	2396S	2396P	D-1396
			SL-2396-I		SL3-90S-I	SL3-90P-I	2396S-I	2396P-I	D-1396-I
		Price	C6		B7	B7	C4	C4	C4
		5	SL-2596		SL5-90S	SL5-90P	2596S	2596P	D-1596
			SL-2596-I		SL5-90S-I	SL5-90P-I	2596S-I	2596P-I	D-1596-I
		Price	C6		C2	C2	C6	C6	C6
		12	MQ-21096		MQ5-90S	MQ5-90P	MQ-21096	MQ-21096	
			MQ-21096-I		MQ5-90S-I	MQ5-90P-I	MQ-21096-I	MQ-21096-I	
					MV5-90S	MV5-90P			
		Price	D3		MQ5-90S-I	MQ5-90P-I	MQ-21096-I	MQ-21096-I	
Welch Allyn	See AAMI								
Protocol	See AAMI								
Micropaq Propaq LT		3	D-1353		D3-90S	D3-90P	2353S	2353P	
			D-1351-I		D3-90S-I	D3-90P-I	2353S-I	2353P-I	
		Price	C7		A8	B3	C1	C1	
		5	D-1553		D5-90S	D5-90P	2553S	2553P	
			D-1553-I		D5-90S-I	D5-90P-I	2553S-I	2553P-I	

PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

		price	C7		B4	B6	C2	C2	
ZOLL									
R SERIES		5	NA		NA	NA	21069S	21069P	
E SERIES							21069S-I	21069P-I	
M SERIES		Price					D3	D3	
M CCT SERIES									

One Piece cables with neonate electrode end for most brands are available upon request. For neonate application, Din Style neonate leads are also the alternatives.

有各品牌新生儿一体化电缆及 D 型电缆+新生儿导联线，使用新生儿心电监护。




1.2. Major types or brands ECG lead wires

主要类型或品牌心电导联线



Types of Major ECG lead wires 主要类型导联线	Lead # 导联	REF# Snap (AHA/IEC) 扣式 (美标/欧标)	OEM REF# OEM 原厂型号	Price code 代格代码	REF# Grabber(AHA/IEC) 钳式 (美标/欧标)	OEM REF# OEM 原厂型号	Price code 代格代码
Din unshielded Individual set leads D 型分离式导联线(非屏蔽)							
	3	D3-90S D3-90S-I		A8	D3-90P D3-90P-I		B3
	5	D5-90S D5-90S-I		B4	D5-90P D5-90P-I		B6
	7	D7-90S/D7-90S-I		B6	D7-90P/D7-90P-I		B7
	10	D10-90S/D10-90S-I		B8	D10-90P/D10-90P-I		C2
Din unshielded individual set (neonate) D 型分离式新生儿导联线 (非屏蔽)							
	3				D3-90NP D3-90NP-I (Mini Grabber)		B4
	5				D5-90NP D5-90NP-I (Mini Grabber)		B6

LL style(shielded) LL 型分离式导联线(屏蔽)							
	3	LL3-90S LL3-90S-I		B7	LL3-90P LL3-90P-I		B7
	5	LL5-90S LL5-90S-I		C2	LL5-90P LL5-90P-I		C2
Mindray/Datascope lead 心电导联线							
	3	MR3-90S MR3-90S-I		C1	MR3-90P MR3-90P-I		C1
	5	MR5-90S MR5-90S-I	012-00-1503-02	C2	MR5-90P MR5-90P-I		C2
	3	DT3-90S DT3-90S-I	0012-00-1503-04/5/6 45cm,61cm,100cm	C2	DT3-90P DT3-90P-I		C2
	5	DT5-90S DT5-90S-I	0012-00-1503-01/2/3 45cm,61cm,100cm	C2	DT5-90P DT5-90P-I	0012-00-1503-01/2/3 45cm,61cm,100cm	C2

Datex-GE Integrated leads (DX300 series) 集成式导联线 (DX300 系列)							
	3	DX3-90S DX3-90S-I	545327	B7	DX3-90P DX3-90P-I	545317 545315	B7
	5	DX5-90S DX5-90S-I	545328	C2	DX5-90P DX5-90P-I	545318 545316	C2
Datex Individual set leads (DX300 series) 分离式导联线 (DX300 系列)							
	3	DXB3-90S DXB3-90S-I	545327	B8	DXB3-90P DXB3-90P-I	545317 545315	B8
	5	DXB5-90S DXB5-90S-I	545328	C4	DXB5-90P DXB5-90P-I	545318 545316	C4
GE Marquette/Multilink-integrated leads 集成式导联线							
	3	MQ3-90S MQ3-90S-I	411203-001,411203-002(130cm) 411203-003,411203-004(130cm)	C2	MQ3-90P MQ3-90P-I	412682-001,412682-002(130cm)) 412682-003,412682-004(130cm))	C2
	5	MQ5-90S MQ5-90S-I	411202-001,411202-002(130cm) 411200-003,411200-004	C2	MQ5-90P MQ5-90P-I	412681-001,412681-002(130cm))	C2

PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

		(Gray)			(Gray)	414556-003,414556-004	
	5	MV5-90S	416468-001/416468-002	C2	MV5-90P	416467-001/416467-002	C2
	V2-V6	MV5-90S-I	416468-003/416468-004	C2	MV5-90P-I	416467-003/416467-004	C2
	C2-C6	(Brown)			(Brown)		
GE Multilink-individual set leads 分离式导联线							
	5	MQB5-90S MQB5-90S-I (Gray)	411200-001,411200-002(130cm) 411200-003,411200-004	C4	MQB5-90P MQB5-90P-I (Gray)	414556-001,414556-002(130cm)) 414556-003,414556-004	C4
	5 V2-V6 C2-C6	MVB5-90S MVB5-90S-I (Brown)	416468-001(130cm),416468-003 416468-004,416468-002	C4	MVB5-90P MVB5-90P-I (Brown)	416467-001(130cm),416467-003 416467-004,416467-002	C4
Nihon Kohden Integrated leads 光电心电导联线							
	3	NKA3-90S NKA3-90S-I	019P	B7	NKA3-90P NKA3-90P-I		B7
	5	NKA5-90S NKA5-90S-I		C2	NKA5-90P NKA5-90P-I		C2
 旧的图片	2	NKB2-90S NKB2-90S-I	BR-912PA BR-912P	NA	NKB2-90P NKB2-90P-I	BR-902PA BR-902P	NA
	3	NKB3-90S NKB3-90S-I	BR-913PA BR-913P	B7	NKB3-90P NKB3-90P-I	BR-903PA BR-903P	B7
	6	NKB6-90S NKB6-90S-I	BR-916PA BR-916P	C3	NKB6-90P NKB6-90P-I	BR-906PA BR-906P	C3

Philips/HP Integrated leads 集成式导联线							
	3	AA3-90S AA3-90S-I	M1605A M1615A	B7	AA3-90P AA3-90P-I	M1603A M1613A	B7
	5	AA5-90S AA5-90S-I	M1615A M1635A	C2	AA5-90P AA5-90P-I	M1623A M1633A	C2
Philips/HP Individual set leads 分离式导联线							
	3	AAB3-90S AAB3-90S-I	M1605A M1615A	B8	AAB3-90P AAB3-90P-I	M1603A M1613A	B8
	5	AAB5-90S AAB5-90S-I	M1615A M1635A	C4	AAB5-90P AAB5-90P-I	M1623A M1633A	C4
Philips IntelliVue Integrated leads 心电导联线							
	3	HP3-90S HP3-90S-I (Gray)	M1673A M1674A	C2	HP3-90P HP3-90P-I (Gray)	M1671A M1672A	C2

PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

	5	HPA5-90S HPA5-90S-I (Blue)	M1644A M1645A	C4	HPA5-90P HPA5-90P-I (Blue)	M1968A M1971A	C4
	5 V2-V6 C2-C6	HPB5-90S HPB5-90S-I (red)	M1602A M1604A	C4	HPB5-90P HPB5-90P-I (red)	M1976A M1978A	C4
Phyiso Control/Meditronic integrated leads							
	4	MDA4-120S MDA4-120S-I		C2			
	6	MDB6-90S MDB6-90S-I		C6			
Siemens/VS,Drager Integrated leads 集成式导联线							
	3	SM3-90S SM3-90S-I		B7	SM3-90P SM3-90P-I		B7
	5	SM5-90S SM5-90S-I		C2	SM5-90P SM5-90P-I		C2
	3 (5/3)	SMB3-90S SMB3-90S-I		C1	SMB3-90P SMB3-90P-I		C1

Siemens/VS,Drager Individual set leads 分离式导联线							
	3	SMD3-90S SMD3-90S-I		B8	SMD3-90P SMD3-90P-I		B8
	5	SMD5-90S SMD5-90S-I		C4	SMD5-90P SMD5-90P-I		C4
	3 (5/3)	SMF3-90S SMF3-90S-I		C1	SMF3-90P SMF3-90P-I		C1
Drager Integrated leads 集成式导联线							
	3	DG3-90S DG3-90S-I		C5	DG3-90P DG3-90P-I		C5
	5	DG5-90S DG5-90S-I		C5	DG5-90P DG5-90P-I		C5

PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

	6	DG6-90S DG6-90S-I		C5	DG6-90P DG6-90P-I		C5
Spacelabs Integrated leads 集成式导联线							
	3	SL3-90S	700-0007-00,60cm	B7	SL3-90P	700-0006-00,60cm	B7
		SL3-90S-I	700-0007-02,100cm		SL3-90P-I	700-0006-02,100cm	
			700-0007-01,60cm			700-0006-01,60cm	
			700-0007-03,100cm			700-0006-03,100cm	
	5	SL5-90S	700-0007-08,60cm	C2	SL5-90P	700-0006-08,60cm	C2
		SL5-90S-I	700-0007-10,100cm		SL5-90P-I	700-0006-10,100cm	
			700-0007-09,60cm			700-0006-09,60cm	
			700-0007-11,100cm			700-0006-11,100cm	
Spacelabs Individual set leads 分离式导联线							
	3	SLB3-90S	700-0007-00,700-0007-02	B8	SLB3-90P	700-0006-00,700-0006-02	B8
		SLB3-90S-I	700-0007-01,700-0007-03		SLB3-90P-I	700-0006-01,700-0006-03	
	4	SLB4-90S	012-0499-06,60cm	C3	SLB4-90P	700-0006-06,100cm	C3
		SLB4-90S-I			SLB4-90P-I	700-0006-07,100cm	
	5	SLB5-90S	700-0007-08,700-0007-10	C4	SLB5-90P	700-0006-08,700-0006-10	C4
		SLB5-90S-I	700-0007-09,700-0007-11		SLB5-90P-I	700-0006-09,700-0006-11	
	6	SLB6-90S		C5	SLB6-90P	700-0006-14,100cm	C5
		SLB6-90S-I			SLB6-90P-I	700-0006-15,100cm	

2. Telemetry ECG Cables 心电遥测电缆



Compatibility 兼容性

Major brands ECG transmitter & Holter system
主要品牌遥测及动态心电系统。



<https://endo.id>

Specification 规格

OEM & Models OEM 原厂设备	Photo 图片	Lead# 导联数	Unimed REF# 尤迈型号	OEM PN 原厂型号	Price code 价格代码
GE Healthcare Apexpro FH Telemetry Patientnet wireless telemetry systems		3	AP3-90S	394111-002	C6
			AP3-90S-I		C6
		5	AP5-90S	394111-006	C7
			AP5-90S-I		C7
		6	AP6-90S	394111-008	C8
			AP6-90S-I		C8
GE Healthcare ApexPro Telemetry		6	MQB6-90S	421930-001	C6
			MQB6-90S-I		
			MQB6-90P	421932-001	C6
		6	MQB6-120S	421930-002	C6
			MQB6-120S-I		

Cable material 线材	TPU Jacket, Light gray TPU 外壳, 灰白色	Connector & Yoke 接头和分线器	Nylon plug & Golden plated pins 尼龙插头和镀金插针
Lead wire length 导联线长度	0.9m		

PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

Mindray/TMS-6012,TEL-100  Mennen/		5	MQB6-120P	421930-002	C6
			TM5-90S		C7
			TM5-90S-I		C7
Philips/HP Medical Legacy M2600A/ M2601A, M2636A Transmitters		3	AT3-90S	M2590A	C7
			AT3-90P	M2591A	
			AT3-90S-I	M2594A	C7
			AT3-90P-I	M2595A	
		5	AT5-90S	M2592A	C8
			AT5-90P	M2593A	
			AT5-90S-I	M2596A	C8
			AT5-90P-I	M2597A	
Telemetry 3/5-lead combiner (U.S. only 仅美国区域)	1 pcs/package	3	Y05-T202	M2598A/M2599A	B8
		5		M2599A	
Philips Medical IntelliVue Telemetry M4841A M2601B		3	HT3-90S	989803151991	C2
			HT3-90P	989803173141	
			HT3-90S-I	989803152001	C2
			HT3-90S-I	989803151981	
		5	HT5-90S	989803173151	C4
			HT5-90P	989803152071	

PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

			HT5-90S-I HT5-90P-I	989803152081 989803152061	C4
Shower Shield (防水罩)	10 pcs/package	3	Y05-T200	989803153031	C3
		5	Y05-T087	989803153041	C3
		6	Y05-T201	989803153051	C4
Single/double ECG alignment guide protective insert, untethered	10 pcs/package	1	Y05-T088	989803140401	A2
		2	Y05-T089	98980314042	A3
MGY 美高仪/MGY-M1 藍港/Visual Stethoscope 3.5mm Jack plug		3	H353-90S H353-90P		B8
			H353-90S-I H353-90P-I		B8
Creative 科瑞康 PC-80 Easy ECG Monitor 2.5mm jack plug		3	H253-90S H253-90P		B8
			H253-90S-I H253-90P-I		B8


3. Disposable ECG lead wires 一次性心电导联线

- Patient Safety
- 安全

- Patient Comfort
- 舒适

- Electrical Performance
- 电气性能




- Ease of Use
- 易使用








Single patient use (disposable) from patient to monitor Reduces the risk of cross contamination from reusable cables and leads, improving patient care and cost control
连接病人与监护仪的一次性心电导联线（供单个病人使用），减少了重复性电缆及导联线在病人间交叉感染的风险，改善了护理环境并有节省成本。

- Full compatibility of the lead wires enable to be used with original ECG trunk Cables or devices, no extra adapter required. 全面兼容原装主电缆或仪器，无需其它转接线
- Built-in ribbon cable for patient comfort and ease of use 置入式带状电缆，舒适易用
- Easy to read labeling facilities easy connection and lead placement 标签清晰易读、连接放置简单
- Unique light blue color indicated the single use 独特的浅蓝色插头设计





3.1 Disposable ECG lead wires 一次性监护分体心电导联线

Descriptions 描述	Probe Connector 导联线接头端	Lead # 导联	REF# Snap(AHA/IEC) 扣式 (美标/欧标)	REF#Grabber(AHA/IEC) 钳式 (美标/欧标)	Packaging 包装	Price code 价格代码
Datascope						
Disposable 3-lead shielded Ribbon cable, 90cm 一次性 3 导联心电导联线, 可分离排线, 屏蔽, 90 厘米		3	DT3-90DS DT3-90DS-I	DT3-90DP DT3-90DP-I	10	B2
Disposable 5-lead shielded Ribbon cable, 90cm 一次性 5 导联心电导联线, 可分离排线, 屏蔽, 90 厘米		5	DT5-90DS DT5-90DS-I	DT5-90DP DT5-90DP-I	10	B2
Drager /Siemens						
Disposable 3-lead shielded Ribbon cable, 90cm 一次性 3 导联心电导联线, 可分离排线, 屏蔽, 90 厘米		3	DG3-90DS DG3-90DS-I	DG3-90DP DG3-90DP-I	10	B3




PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

Disposable 5-lead shielded Ribbon cable, 90cm 一次性 5 导联心电图导联线，可分离排线，屏蔽，90 厘米		5	DG5-90DS DG5-90DS-I	DG5-90DP DG5-90DP-I	10	B3
Disposable 6-lead shielded Ribbon cable, 90cm 一次性 6 导联，可分离排线，屏蔽，90cm		6	DG6-90DS DG6-90DS-I	DG6-90DP DG6-90DP-I	10	B5
Disposable 3-lead shielded Ribbon cable, 90cm 一次性 3 导联心电图导联线，可分离排线，屏蔽，90 厘米		3	SM3-90DS SM3-90DS-I	SM3-90DP SM3-90DP-I	10	B2
Disposable 5-lead shielded Ribbon cable, 90cm 一次性 5 导联心电图导联线，可分离排线，屏蔽，90 厘米		5	SM5-90DS SM5-90DS-I	SM5-90DP SM5-90DP-I	10	B2
GE Multi-link						
Disposable 3-lead shielded Ribbon cable, 90cm 一次性 3 导联，可分离排线，屏蔽，90 厘米		3	MQ3-90DS MQ3-90DS-I	MQ3-90DS MQ3-90DS-I	10	B2
Disposable 5-lead shielded Ribbon cable, 90cm 一次性 5 导联，可分离排线，屏蔽，90cm		5	MQ5-90DS MQ5-90DS-I	MQ5-90DP MQ5-90DP-I	10	B2




PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

Disposable 5-lead shielded Ribbon cable, V or C lead ,90cm 一次性 5 导联，可分离排线，V 或 C 导联，屏蔽，90cm		5	MV5-90DS MV5-90DS-I	MV5-90DP MV5-90DP-I	10	B2
Disposable 6-lead shielded Ribbon cable for Apex, 90cm 一次性 6 导联，可分离排线，屏蔽，用于 APEX，90cm		6	MQB6-90DS MQB6-90DS-I	MQB6-90DP MQB6-90DP-I	10	B5
Mindray						
Disposable 3-lead shielded Ribbon cable, 90cm 一次性 3 导联，可分离排线，屏蔽，90 厘米		3	MR3-90DS MR3-90DS-I	MR3-90DP MR3-90DP-I	10	B2
Disposable 5-lead shielded Ribbon cable, 90cm 一次性 5 导联导联线，可分离排线，屏蔽，90 厘米		5	MR5-90DS MR5-90DS-I	MR5-90DS MR5-90DS-I	10	B2
Nihon Kohden						
Disposable 3-lead shielded Ribbon cable, 90cm 一次性 3 导联心电导联线，可分离排线，屏蔽，90 厘米		3	NKB3-90DS NKB3-90DS-I	NKB3-90DP NKB3-90DP-I	10	B2
Disposable 6-lead shielded Ribbon cable, 90cm 一次性 6 导联，可分离排线，屏蔽，90cm		6	NKB6-90DS NKB6-90DS-I	NKB6-90DP NKB6-90DP-I	10	B5
Philips IntelliVue						


PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

Disposable 3-lead shielded Ribbon cable, 90cm 一次性 3 导联心电图导联线，可分离排线，屏蔽，90 厘米		3	HP3-90DS HP3-90DS-I	HP3-90DP HP3-90DP-I	10	B3
Disposable 5-lead shielded Ribbon cable, 90cm, 一次性 5 导联心电图导联线，可分离排线，屏蔽，90cm		5	HPA5-90DS HPA5-90DS-I	HPA5-90DS-I HPA5-90DP-I	10	B3
Disposable 5-lead shielded Ribbon V or C-lead cable, 90cm, 一次性 5 导联心电图导联线，可分离排线，屏蔽，V 或 C-导，90cm		5	HPB5-90DS HPB5-90DS-I	HPB5-90DP HPB5-90DP-I	10	B3
Disposable 3-lead shielded Ribbon cable, 90cm 一次性 3 导联心电图导联线，可分离排线，屏蔽，90 厘米		3	HT3-90DS HT3-90DS-I	HT3-90DP HT3-90DP-I	10	B5
Disposable shielded Ribbon cable with shower shielded, 90cm 一次性 5 导联心电图导联线，可分离排线，屏蔽，带防水罩，90cm		5	HT5-90DS HT5-90DS-I	HT5-90DP HT5-90DP-I	10	B5
Philips Twin Pin(AA style)						

PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

Disposable 3-lead shielded Ribbon cable, 90cm 一次性 3 导联，可分离排线，屏蔽，90cm		3	AAB3-90DS AAB3-90DS-I	AAB3-90DP AAB3-90DP-I	10	B3
Disposable 5-lead shielded Ribbon cable, 90cm 一次性 5 导联，可分离排线，屏蔽，90cm		5	AAB5-90DS AAB5-90DS-I	AAB5-90DP AAB5-90DP-I	10	B3
Philips Veridia Telemetry						
Disposable 3-lead shielded Ribbon cable, 90cm with combiner 一次性 3 导联导联线，带线框，可分离排线，屏蔽，90cm		3	AT3-90DS AT3-90DS-I	AT3-90DP AT3-90DP-I	10	B5
Disposable 5-lead shielded Ribbon cable, 90cm with combiner 一次性 3 导联导联线，带线框，可分离排线，屏蔽，90cm		5	AT5-90DS AT5-90DS-I	AT5-90DP AT5-90DP-I	10	B5
Spacelabs						
Disposable 3-lead shielded Ribbon cable, 90cm 一次性 3 导联，可分离排线，屏蔽，90 厘米	NA	3	SLB3-90DS SLB3-90DS-I	SLB3-90DP SLB3-90DP-I	10	B2

PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

Disposable 5-lead shielded Ribbon cable, 90cm 一次性 5 导联导联线，可分离排线，屏蔽，90 厘米		5	SL5-90DS SL5-90DS-I	SL5-90DP SL5-90DP-I	10	B2
---	--	---	--------------------------------------	--------------------------------------	----	----


3.2 Disposable ECG One piece lead wires 一次性诊断心电导联线

OEM & Models 原厂及型号	EKG machine Plug 机器接头端	Lead # 导联数	REF# Snap(AHA/IEC) 扣式（美标/欧标）	REF# Grabber(AHA/IEC) 钳式（美标/欧标）	Packaging 包装	Unit code 单根价格
AAMI						
		3	2340DS 2340DS-I	2340DP 2385DP-I	8	B6.5
		5	2540DS 2540DS-I	2540DP 2540DP-I	8	B7
Philips						
Merlin 78352C/78354C/ 78834CM1001A/		3	2385DS 2385DS-I	2385DP 2385DP-I	8	B7
		5	2585DS 2585DS-I	2585DP 2585DP-I	8	C1.3
GE Multi-link						
Eagle,Solar,Dash monitor MAC-Lab Cath Lab		3	2386DS 2386DS-I	2386DP 2386DP-I	8	B6.5

PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

Datex-Ohmeda M-Prestn & “E” modules		5	2586DS 2586DS-I	2586DP 2586DP-I	8	B7
Nihon Kohden						
BSM-2300 lifescope I BSM-4100 lifescope P BSM-5100		3	2309DS 2309DS-I	2309DP 2309DP-I	8	B7
Mindray / Datascope						
PM5000 PM6000		3	2312DS 2312DS-I	2312DP 2312DP-I	8	B6.5
		5	2512DS 2512DS-I	2512DP 2512DP-I	8	B7
Datascope -Mindray V12, V21		3	2326DS 2326DS-I	2326DP 2326DP-I	8	B7
		5	2526DS 2512DS-I	2526DP 2526DP-I	8	C1.3
Spacelabs						
Ultraview 90492, 91220, 90496, 91496		3	2396DS 2396DS-I	2396DP 2396DP-I	8	B7
		5	2596DS 2596DS-I	2596DP 2596DP-I	8	C1.3

4. ECG Lead Adapter 心电导联转换头



Compatibility 兼容:





Major brands ECG monitoring leads system

The unique design of the ECG lead wires adaptors enable the hospitals to use the **one type lead wires** to fit major brands patient monitoring system running in one hospital by attaching the proper ECG adaptor to the trunk cable. No confusion on lead wires any more!

心电导联线转换器的独特设计，让医院实现了用一种导联线适配各种主要品牌医用监护系统的可能。只须让合适的转换器连上主电缆，一切关于导连线的选用困扰都将不复存在！



Specification 产品规格:

<https://endo.id>

OEM Manufacturer OEM 制造商	3-lead REF# 3-导联型号	Trunk Cable end 主电缆端接头	Lead wire end 导联线端接头	Price 价格	5-lead REF# 5-导联型号	Trunk Cable end 主电缆端接头	Lead wire end 导联线端接头	Price code 价格代码
Convert To Din style lead wires D 型导联线转接头								
Colin	AD-CL3-D3			B5	AD-CL5-D5			B5
GE Marquette	AD-MQ5-D3			B5	AD-MQ5-D5			B5
GE Datex-Ohmeda	AD-DX3-D3			B5	AD-DX5-D5			B5
Nihon Kohden	AD-NKA3-D3			B5	AD-NKA5-D5			B5
Nihon Kohden	AD-NKB3-D3			B5				
Philips medical	AD-AA3-D3			B5				
Philips medical	AD-HP3-D3			B8				
Siemens/Drager	AD-SM3-D3			B5	AD-SM5-D5			B5
Spacelabs	AD-SL3-D3			B5	AD-SL5-D5			B5




Cable 线材		Connector & Yoke 接头和分线器	Nylon plug & Golden plated pins 尼龙插头和镀金插针
Link 绑带			
http			





OEM Manufacturer OEM 制造商	3-lead REF# 3-导联型号	Trunk Cable end 主电缆端接头	Lead wire end 导联线端接头	Price 价格	5-lead REF# 5-导联型号	Trunk Cable end 主电缆端接头	Lead wire end 导联线端接头	Price code 价格代码
Convert To VS/Siemens lead	VS/Siemens 导联线转接头							
Colin	AD-CL3-SM3			B5	AD-CL5-SM5			B5
GE Marquette	AD-MQ5-SM3			B5	AD-MQ5-SM5			B5
GE Datex-Ohmeda	AD-DX3-SM3			B5	AD-DX5-SM5			B5
Nihon Kohden	AD-NKA3-SM3			B5	AD-NKA5-SM5			B5
Nihon Kohden	AD-NKB3-SM3			B5				
Philips medical	AD-AA3-SM3			B5				
Philips medical	AD-HP3-SM3			B8				
Spacelabs	AD-SL3-SM3			B5	AD-SL5-SM5			B5

PT. ENDO Indonesia -- Importer, distributor, and manufacturer of



OEM Manufacturer OEM 制造商	3-lead REF# 3-导联型号	Trunk Cable end 主电缆端接头	Lead wire end 导联线端接头	Price 价格	5-lead REF# 5-导联型号	Trunk Cable end 主电缆端接头	Lead wire end 导联线端接头	Price code 价格代码
Convert To Philips/HP style lead		Philips/HP 型转换头						
Colin	AD-CL3-HP3			B6	AD-CL5-HP5			B7
GE Marquette	AD-MQ5-HP3			B6	AD-MQ5-HP5			B7
GE Datex-Ohmeda	AD-DX3-HP3			B6	AD-DX5-HP5			B7
Nihon Kohden	AD-NKA3-HP3			B6	AD-NKA5-HP5			B7
Nihon Kohden	AD-NKB3-HP3			B6				
Philips medical	AD-AA3-HP3			B6				
Siemens/Drager	AD-SM3-HP3			B6	AD-SM5-HP5			B7
Spacelabs	AD-SL3-HP3			B6	AD-SL5-HP5			B7






OEM Manufacturer OEM 制造商	3-lead REF# 3-导联型号	Trunk Cable end 主电缆端接头	Lead wire end 导联线端接头	Price 价格	5-lead REF# 5-导联型号	Trunk Cable end 主电缆端接头	Lead wire end 导联线端接头	Price code 价格代码
Convert To GE Multi-Link lead	GE Multi-Link 转换头							
Colin					AD-CL5-MQ5			B5
GE Datex-Ohmeda					AD-DX5-MQ5			B5
Nihon Kohden					AD-NKA5-MQ5			B5
Siemens/Drager					AD-SM5-MQ5			B5
Spacelabs					AD-SL5-MQ5			B5






Application Figure 应用图片



PT ENDO Indonesia is a distributor, and manufacturer of medical equipment.



OEM Manufacturer OEM 制造商	3-lead REF# 3-导联型号	Trunk Cable end 主电缆端接头	Lead wire end 导联线端接头	Price 价格	5-lead REF# 5-导联型号	Trunk Cable end 主电缆端接头	Lead wire end 导联线端接头	Price code 价格代码
Convert To Datex 300 style lead								
Colin	AD-CL3-DX3			B5	AD-CL5-DX5			B5
GE Marquette					AD-MQ5-DX5			B5
Nihon Kohden	AD-NKA3-DX3			B5	AD-NKA5-DX5			B5
Nihon Kohden	AD-NKB3-DX3			B5				
Philips medical	AD-AA3-DX3			B5				
Philips medical	AD-HP3-DX3			B6				
Siemens/Drager	AD-SM3-DX3			B5	AD-SM5-DX5			B5
Spacelabs	AD-SL3-DX3			B5	AD-SL5-DX5			B5

OEM Manufacturer OEM 制造商	3-lead REF# 3-导联型号	Trunk Cable end 主电缆端接头	Lead wire end 导联线端接头	Price 价格	5-lead REF# 5-导联型号	Trunk Cable end 主电缆端接头	Lead wire end 导联线端接头	Price code 价格代码
Convert To Spacelabs style lead								
Colin	AD-CL3-SL3			B5	AD-CL5-SL5			B5
GE Marquette	AD-MQ5- SL3			B5	AD-MQ5-SL5			B5
GE Datex-Ohmeda	AD-DX3- SL3			B5	AD-DX5-SL5			B5
Nihon Kohden	AD-NKA3- SL3			B5	AD-NKA5-SL5			B5
Nihon Kohden	AD-NKB3- SL3			B5				
Philips medical	AD-AA3- SL3			B5				
Philips medical	AD-HP3- SL3			B6				
Siemens/Drager	AD-SM3- SL3			B5	AD-SM5-SL5			B5

5.Diagnosis ECG Cables & Accessories 诊断用心电电缆及附件

Specification 产品规格

Compatibility

兼容性

Major brands ECG diagnosis system.

主要品牌的心电图机，运动平板等心电诊断设备。



Product name 产品类型	Photo 图片	Part # Example 型号#例举	R
-------------------	----------	----------------------	---







Cable material 线材	TPU Jacket, Light gray TPU 外被, 浅灰	Connector 接头	Nylon plug & Golden plated pins 尼龙插头和镀金插针
Main Cable length 主电缆长度	2.2m	Lead wire length 导联线长度	1.2m(limb lead);0.8m(Chest lead) 1.2m (四肢); 0.8m(胸腔)

<https://endo.id>

One piece Cable 一体化电缆		E10R-HP-B/I (Example 例如)	E: EKG diagnostic cable or lead wire E:心电图机电缆或导联线 10: 10-lead(12=12 lead) 10: 10 导(12=12 导) R: built in resistors for protection(default: no resistors) R: 带电阻(缺省: 无电阻) HP: Philips/HP,factory/model code HP: Philips /HP,原厂/机型 1:Connector configure; 1:接头 类型; B: banana end B: 香蕉头 (N= needle, P=grabber, S=snap) (N= 针式, P=钳式, S=扣式) I: IEC code (Default is AHA code) I: 欧标 (缺省:美标)
Trunk Cable 主电缆		E10R-MQ-I (Example 例如)	E: EKG diagnostic cable or lead wire E:心电图机电缆或导联线 10: 10-lead(12=12 lead) 10: 10 导(12=12 导) R: built in resistors for protection(default: no resistors) R: 带电阻(缺省: 无电阻) MQ: GE Marquette, factory/model code MQ: GE Marquette, 原厂/机型 1:Connector configure (default: only one connector available) 1: 接头类型 (缺省: 仅一种接头, 不可选) I: IEC code (Default is AHA code) I: 欧标 (缺省:美标)
Lead wire set 导联线		MQ10-LB-I (Example 例如)	MQ: Marquette, factory/model code MQ: GE Marquette,原厂/机型 10: 10: 10-lead(12=12 lead) 10: 10 导(12=12 导) L: lead wires L: 导联线 B: banana end B: 香蕉头 (N= needle, P=grabber, S=snap) (N= 针式, P=钳式, S=扣式) I: IEC code (Default is AHA code) I: 欧标 (缺省:美标)

PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia







All one piece cable or lead wires are available with banana end, needle end, grabber end and snap end 供应各式一体化电缆和导联线 (香蕉头 / 针式 / 钳式 / 扣式)

Code 代码	AHA 美标	IEC 欧标	IH 替代欧标
4mm Banana end 香蕉头 (B)			
3mm Needle end 针式 (N)			

Code 代码	AHA 美标	IEC 欧标	IH 替代欧标
---------	--------	--------	---------

<https://endo.id>




PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

Grabber End 钳式 (P)			
Snap End 扣式 (S)			





5.1. One piece EKG fixed cable 一体化心电电缆

<https://endo.id>

PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

OEM& Models 原厂设备	EKG machine Plug 机器端插头	Lead# 导联数	With 4mm banana 一体化香蕉头 4mm	With 3mm needle 一体化针头 3mm	With grabber end 一体化钳式	With Snap end 一体化扣式	Price code 价格代码
Burdick/Cardiac Science							
EK-10 and Elite, E350, E350i, E550, E560, Eclipse™ LE/LE.		10	E10-BK1-B E10R-BK1-B (012-0700-00)	E10-BK1-N E10R-BK1-N	E10-BK1-P E10R-BK1-P	E10-BK1-S E10R-BK1-S	D3
		10	E10-BK1-B/I E10R-BK1-B/I	E10-BK1-N/I E10R-BK1-N/I	E10-BK1-P/I E10R-BK1-P/I	E10-BK1-S/I E10R-BK1-S/I	D3
Burdick Quinton Eclipse 4,8,400,800, 850,850i,plus, Eclipse LE II and ATRIA PLUS Quinton?		10	E10-BK2-B E10R-BK2-B (007704)	E10-BK2-N E10R-BK2-N	E10-BK2-P E10R-BK2-P	E10-BK2-S E10R-BK2-S	D3
		10	E10-BK2-B/I E10R-BK2-B/I	E10-BK2-N/I E10R-BK2-N/I	E10-BK2-P/I E10R-BK2-P/I	E10-BK2-S/I E10R-BK2-S/I	D3
Fukuda Denshi							
FX101, FX102F, FX328U, FCP2201		10	E10-FD1-B E10R-FD1-B	E10-FD1-N E10R-FD1-N	E10-FD1-P E10R-FD1-P	E10-FD1-S E10R-FD1-S	D3
		10	E10-FD1-B/I E10R-FD1-B/I	E10-FD1-N/I E10R-FD1-N/I	E10-FD1-P/I E10R-FD1-P/I	E10-FD1-S/I E10R-FD1-S/I	D3
Fukuda ME							
KP-500、KP-500D		10	E10-FD2-B E10R-FD2-B	E10-FD2-N E10R-FD2-N	E10-FD2-P E10R-FD2-P	E10-FD2-S E10R-FD2-S	D3
		10	E10-FD2-B/I E10R-FD2-B/I	E10-FD2-N/I E10R-FD2-N/I	E10-FD2-P/I E10R-FD2-P/I	E10-FD2-S/I E10R-FD2-S/I	D3




PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

OEM& Models 原厂设备	EKG machine Plug 机器端插头	Lead# 导联数	With 4mm banana 一体化香蕉头 4mm	With 3mm needle 一体化针头 3mm	With grabber end 一体化钳式	With Snap end 一体化扣式	Price code 价格代码
GE-Marquette							
MAC 500, MAC 1100, MAC 1200		10	E10-MQ-B E10R-MQ-B	E10-MQ-N E10R-MQ-N	E10-MQ-P E10R-MQ-P	E10-MQ-S E10R-MQ-S	D3
		10	E10-MQ-B/I E10R-MQ-B/I	E10-MQ-N/I E10R-MQ-N/I	E10-MQ-P/I E10R-MQ-P/I	E10-MQ-S/I E10R-MQ-S/I	D3
Hellige / Siemens Hörmann / Bosch							
		10	E10-HL-B E10R-HL-B	E10-HL-N E10R-HL-N	E10-HL-P E10R-HL-P	E10-HL-S E10R-HL-S	D3
		10	E10-HL-B/I E10R-HL-B/I	E10-HL-N/I E10R-HL-N/I	E10-HL-P/I E10R-HL-P/I	E10-HL-S/I E10R-HL-S/I	D3
Kenz							
103,106		10	E10-KZ1-B E10R-KZ1-B	E10-KZ1-N E10R-KZ1-B	E10-KZ1-P E10R-KZ1-P	E10-KZ1-S E10R-KZ1-S	D3
CARDIOLINE (Italy)		10	E10-KZ1-B/I E10R-KZ1-B/I	E10-KZ1-N/I E10R-KZ1-N/I	E10-KZ1-P/I E10R-KZ1-P/I	E10-KZ1-S/I E10R-KZ1-S/I	D3
ECG 108/110 Cardico 302/601/1210/1211		10	E10-KZ2-B E10R-KZ2-B	E10-KZ-N E10R-KZ-B	E10-KZ2-P E10R-KZ2-P	E10-KZ2-S E10R-KZ2-S	D3

<https://endo.id>



PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

		10	E10-KZ2-B/I E10R-KZ2-B/I	E10-KZ2-N/I E10R-KZ2-N/I	E10-KZ2-P/I E10R-KZ2-P/I	E10-KZ2-S/I E10R-KZ2-S/I	D3
OEM& Models 原厂设备	EKG machine Plug 机器端插头	Lead# 导联数	With 4mm banana 一体化香蕉头 4mm	With 3mm needle 一体化针头 3mm	With grabber end 一体化钳式	With Snap end 一体化扣式	Price code 价格代码
Mortara							
ELI-200、ELI-50		10	E10-MT-B E10R-MT-B	E10-MT-N E10R-MT-N	E10-MT-P E10R-MT-P	E10-MT-S E10R-MT-S	D3
		10	E10-MT-B/I E10R-MT-B/I	E10-MT-N/I E10R-MT-N/I	E10-MT-P/I E10R-MT-P/I	E10-MT-S/I E10R-MT-S/I	D3
Nihon Kohden							
CaduoFax 6151,6851,6353 CadioFax 8110		10	E10-NK1-B E10R-NK1-B	E10-NK1-N E10R-NK1-N	E10-NK1-P E10R-NK1-P	E10-NK1-S E10R-NK1-S	D3
		10	E10-NK1-B/I E10R-NK1-B/I	E10-NK1-N/I E10R-NK1-N/I	E10-NK1-P/I E10R-NK1-P/I	E10-NK1-S/I E10R-NK1-S/I	D3
CaduoFax Q ECG-9110K ECG-9130K Doing Juang Cardioline (New)		10	E10-NK2-B E10R-NK2-B	E10-NK2-N E10R-NK2-N (BA-901D)	E10-NK2-P E10R-NK2-P	E10-NK2-S E10R-NK2-S	D3
		10	E10-NK2-B/I E10R-NK2-B/I (BJ-902D)	E10-NK2-N/I E10R-NK2-N/I (BJ-901D)	E10-NK2-P/I E10R-NK2-P/I	E10-NK2-S/I E10R-NK2-S/I	D3
Philips/HP							
M1770A,M1771A,M1772A Pagewrite100/200/ 200i		10	E10-HP-B E10R-HP-B (M3702C)	E10-HP-N E10R-HP-N	E10-HP-P E10R-HP-P	E10-HP-S E10R-HP-S	D3



PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

		10	E10-HP-B/I E10R-HP-B/I (M3703A/M3703C)	E10-HP-N/I E10R-HP-N/I	E10-HP-P/I E10R-HP-P/I	E10-HP-S/I E10R-HP-S/I	D3
OEM& Models 原厂设备	EKG machine Plug 机器端插头	Lead# 导联数	With 4mm banana 一体化香蕉头 4mm	With 3mm needle 一体化针头 3mm	With grabber end 一体化钳式	With Snap end 一体化扣式	Price code 价格代码
Schiller							
AT1, AT2, AT2 Plus, AT4 Resistor 10K Long Screw		10	E10-SH1-B E10R-SH1-B	E10-SH1-N E10R-SH1-N	E10-SH1-P E10R-SH1-P	E10-SH1-S E10R-SH1-S	D3
		10	E10-SH1-B/I E10R-SH1-B/I	E10-SH1-N/I E10R-SH1-N/I	E10-SH1-P/I E10R-SH1-P/I	E10-SH1-S/I E10R-SH1-S/I	D3
AT3, AT6, CS6, AT5, AT10, AT60 Short Screw		10	E10-SH2-B E10R-SH2-B	E10-SH2-N E10R-SH2-N	E10-SH2-P E10R-SH2-P	E10-SH2-S E10R-SH2-S	D3
		10	E10-SH2-B/I E10R-SH2-B/I	E10-SH2-N/I E10R-SH2-N/I	E10-SH2-P/I E10R-SH2-P/I	E10-SH2-S/I E10R-SH2-S/I	D3











5.2. EKG Trunk Cable 心电主电缆

OEM& Models 原厂设备	EKG machine Plug 机器端插头	Lead# 导联数	REF# 型号	Used with lead wire set/GE Multi-link set 与导联线/GE Multi-link 导联线一起使用	Price code 价格代码
Burdick					
EK10-E350		10	E10-BK1 E10R-BK1	10-lead, chest /0.9m;Limb/1.2m MQ10-LB,MQ-LN,MQ-LP,MQ-LS or MQB10-LB,...	D3
		10	E10-BK1-I E10R-BK1-I	10-lead, chest /0.9m;Limb/1.2m MQ10-LB-I,MQ-LN-I,MQ-LP-I,MQ-LS-I or MQB-LB,.... (See the next page of 10-lead wire sets for detail) (详细细节见下页 10-导 导联线)	D3
Burdick Quinton Eclipse 4,8,400,800,		10	E10-BK2 E10R-BK2	Same as above 同上	D3



PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

850,850i,plus, Eclipse LE II and ATRIA PLUS		10	E10-BK2-I E10R-BK2-I	Same as above	D3
Fukuda Denshi					
FX101, FX102F, FX328U, FCP2201		10	E10-FD1 E10R-FD1	Same as above	D3
		10	E10-FD1-I E10R-FD1-I	Same as above	D3
Fukuda ME					
KP-500、KP-500D		10	E10-FD2 E10R-FD2	Same as above	D3
		10	E10-FD2-I E10R-FD2-I	Same as above	D3
OEM& Models 原厂设备	EKG machine Plug 机器端插头	Lead# 导联数	REF# 型号	Used with lead wire set/GE Multi-link set 与导联线/GE Multi-link 导联线一起使用	Price code 价格代码
GE-Marquette					
MAC 500, MAC 1100, MAC 1200		10	E10-MQ E10R-MQ (22341809)	Same as above	D3
		10	E10-MQ-I E10R-MQ-I	Same as above	D3
ELI-200、ELI-50		10	E10-MT E10R-MT	Same as above	D3

PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia


		10	E10-MT-I E10R-MT-I	Same as above	D3
Hellige / Siemens Hörmann / Bosch					
		10	E10-HL E10R-HL	Same as above	D3
		10	E10-HL-I E10R-HL-I	Same as above	D3
Kenz					
103,106		10	E10-KZ1 E10R-KZ1	Same as above	D3
CARDIOLINE (Italy)		10	E10-KZ1-I E10R-KZ1-I	Same as above	D3
ECG 108/110 Cardico 302/601/1210/1211		10	E10-KZ2 E10R-KZ2	Same as above	D3
		10	E10-KZ2-I E10R-KZ2-I	Same as above	D3
OEM& Models 原厂设备	EKG machine Plug 机器端插头	Lead# 导联数	REF# 型号	Used with lead wire set/GE Multi-link set 与导联线/GE Multi-link 导联线一起使用	Price code 价格代码
Nihon Kohden					
CaduoFax 6151,6851,6353 CadioFax 8110		10	E10-NK1 E10R-NK1	Same as above	D3

PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

		10	E10-NK1-I E10R-NK1-I	Same as above	D3
CaduoFax Q ECG-9110K ECG-9130K Doing Juang Cardioline (New)		10	E10-NK2 E10R-NK2	Same as above	D3
		10	E10-NK2-I E10R-NK2-I	Same as above	D3
Philips/HP					
M1770A,M1771A,M1772A Pagewrite100/200/ 200i		10	E10-HP E10R-HP	Same as above	D3
		10	E10-HP-I E10R-HP-I	Same as above	D3
OEM& Models 原厂设备	EKG machine Plug 机器端插头	Lead# 导联数	REF# 型号	Used with lead wire set/GE Multi-link set 与导联线/GE Multi-link 导联线一起使用	Price code 价格代码
Schiller					
AT1, AT2, AT2 Plus,AT4 Resistor 10K		10	E10-SH1 E10R-SH1	Same as above	D3

<https://endo.id>

PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

Long Screw		10	E10-SH1-I E10R-SH1-I	Same as above	D3
AT3,AT6,CS6,AT5, AT10,AT60 Short Screw		10	E10-SH2 E10R-SH2	Same as above	D3
		10	E10-SH2-I E10R-SH2-I	Same as above	D3
GE Healthcare-Resting sys.					
CAM14 coiled patient Cable Used with GE MAC 5000 resting system with CAM14 Acquisition Module		10	ECAM-GE14 (2016560-001) 4.5M+0.65M ECAM-GE14-S (2016560-001) 1.4M+0.3M	Used with CAM 14 patient lead wires on next page. 与下页的 CAM 14 导联线配合使用	D8
CASE 12, CASE 15, MAC PC, MAC 6, MAC 12, MAC 15, CENTRA Systems with AM4 or AM5 Acquisition Module		10	EAM-GE4-S (700044-202)0.25+4.6M EAM-GE4-S (700044-203) 0.25+1.3M	Used with CAM 14patient lead wires on next page 与下页的 CAM 14 导联线配合使用	D8



5.3 EKG Lead wire 心电导联线

<https://endo.id>

PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

GE Multi-Link 10-Lead-integrated leads 集成式导联线	OEM & Models 原厂型号	Electrode end 病人电极端	REF# (AHA) 型号 (美标)	REF# (IEC) 型号 (欧标)	Pieces/Set 件/套	Price code 价格代码
	Banana End 900112-001(AHA) 900112-002(IEC)	4mm banana	MQ10-LB (900112-001)	MQ10-LB-I (900112-002)	10	D1
		3mm needle	MQ10-LN	MQ10-LN-I	10	D1
		Grabber	MQ10-LP	MQ10-LP-I	10	D1
		Snap	MQ10-LS	MQ10-LS-I	10	D1
GE CAM 14 patient lead wire 导联线	OEM & Models 原厂型号	Electrode end 病人电极端	REF# (AHA) 型号 (美标)	REF# (IEC) 型号 (欧标)	Pieces/Set 件/套	Price code 价格代码
		4mm banana	MQB10-LB	MQB10-LB-I	10	D1
		3mm needle	MQB10-LN	MQB10-LN-I	10	D1
		Grabber	MQB10-LP	MQB10-LP-I	10	D1
		Snap	MQB10-LS	MQB10-LS-I	10	D1
	420101-002	CAM/Universal	MQB10-AD	MQB10-AD-I	10	D1
	420101-001 Used with 900178-001/900178-002 or Unimed Adapter as electrodes end (see below)	CAM/Universal	MQB14-AD	MQB14-AD-I	14	D3

<https://endo.id>

PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

	2001925-005(66cm)	CAM/Universal	X-MQB-66DF	66cm,no color coded	10	C8
	2001925-002(74cm)	CAM/Universal	X-MQB-80DF	80cm,no color coded	10	C8
	2001925-006(90cm)	CAM/Universal	X-MQB-90DF	90cm,no color coded	10	C8
	2001925-003(102cm)	CAM/Universal	X-MQB-100DF	100cm,no color coded	10	C8
	2001925-004(130cm)	CAM/Universal	X-MQB-120DF	120cm,no color coded	10	C8
GE AM4/AM5 patient lead wires 导联线	OEM & Models 原厂型号	Electrode end 病人电极端	REF# (AHA) 型号（美标）	REF# (IEC) 型号（欧标）	Pieces/Set 件/套	Price code 价格代码
	?	4mm banana	MQD10-LB	MQD10-LB-I	10	D1
		3mm needle	MQD10-LN	MQD10-LN-I	10	D1
		Grabber	MQD10-LP	MQD10-LP-I	10	D1
		Snap	MQD10-LS	MQD10-LS-I	10	D1
	900177-001 Used with 900178-001/900178-002 or Unimed Adapter as electrodes end (see below)	CAM/Universal	MQD12-AD	MQD12-AD-I	12	D3

PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

	406581-004(90cm)	CAM/Universal	X-MQD-80DF	80cm,no color coded	10	C8
	406581-006(130cm)	CAM/Universal	X-MQD-120DF	120cm,no color coded	10	C8
GE CAM14/AM4/AM5 Electrode Adapter 电极转接头	OEM & Models 原厂型号	Electrode end 病人端电极	REF# (AHA) 型号 (美标)	REF# (IEC) 型号 (欧标)	Pieces/Set 件/套	Price code 价格代码
	900178-001,AHA	4mm banana,10 pcs/set	AD10-DB	AD10-DB-I	10	B8
	900178-101,IEC					
	900179-201,AHA	4mm banana,14 pcs/set	AD14-DB	AD14-DB-I	14	C3
	900178-201,IEC					
	For the above leads set Or GE original CAM14/AM4/AM5	4mm banana,12 pcs/set	AD12-DB	AD12-DB-I	12	C2
Philips Lead wire 导联线	OEM & Models 原厂型号	Electrode end 病人端电极	REF# (AHA) 型号 (美标)	REF# (IEC) 型号 (欧标)	Pieces/Set 件/套	Price code 价格代码
	PageWriter Touch: 860284,	4mm banana	PW10-LB	PW10-LB-I	10	D1
	PageWriter		(989803129161)	(989803129191)		
	Trim Cardiographs, 860286, 860287,	3mm needle	PW10-LN	PW10-LN-I	10	D1
	860288, 860290, 860297	Grabber	PW10-LP	PW10-LP-I	10	D1
		Snap	PW10-LS	PW10-LS-I	10	D1

PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

	989803151631/989803151641	4mm banana	PWC10-LB	PWC10-LB-I	10	D1	
		3mm needle	PWC10-LN	PWC10-LN-I	10	D1	
		Grabber	PWC10-LP	PWC10-LP-I	10	D1	
		Snap	PWC10-LS	PWC10-LS-I	10	D1	
QULNTON Lead wire 导联线							
	Q-710 Stress Test System 031257-002	4mm banana	QT10-LB	QT10-LB-I	10	D1	
		3mm needle	QT10-LN	QT10-LN-I	10	D1	
		Grabber	QT10-LP	QT10-LP-I	10	D1	
		Snap	QT10-LS	QT10-LS-I	10	D1	
Mortara Lead wire 导联线							
	9293-046-60	4mm banana	MT10-LB	MT10-LB-I	10	D1	
		3mm needle	MT10-LN	MT10-LN-I	10	D1	
		Grabber	MT10-LP	MT10-LP-I	10	D1	
		Snap	MT10-LS	MT10-LS-I	10	D1	
Norav (AA style) Lead wire 导联线		OEM & Models 原厂型号	Electrode end 病人端电极	REF# (AHA) 型号（美标）	REF# (IEC) 型号（欧标）	Pieces/Set 件/套	Price code 价格代码
		4mm banana	AA10-LB	AA10-LB-I	10	D1	
		3mm needle	AA10-LN	AA10-LN-I	10	D1	
		Grabber	AA10-LP	AA10-LP-I	10	D1	
		Snap	AA10-LS	AA10-LS-I	10	D1	

<https://endo.id>



PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

Welch Allyn Lead wire 导联线						
	CP100/CP200 401129	4mm banana	WA10-LB	WA10-LB-I	10	D1
		3mm needle	WA10-LN	WA10-LN-I	10	D1
		Grabber	WA10-LP	WA10-LP-I	10	D1
		Snap	WA10-LS	WA10-LS-I	10	D1

5.4 EKG Accessories 附件


ECG Electrodes 心电电极	REF# 型号	Package 包装	Price code 价格代码	REF to Figures 图例
Adult Limb Clamps/IEC, universal 成人四肢夹/欧标, 通用	UU-LC-A	4 pcs/box	B6	A
Adult Limb Clamps/AHA, all clamps in green 成人四肢夹/美标, 绿色	UU-LC-AA	4 pcs/box	B6	E
Pediatric Limb Clamps/IEC, universal 儿童四肢夹/欧标, 通用	UU-LC-P	4 pcs/box	B6	B
Pediatric Limb Clamps/AHA, all clamps in green 儿童四肢夹/美标, 绿色	UU-LC-PA	4 pcs/box	B6	未显示
Adult Limb Clamps, for snap electrode 成人四肢夹, 接扣式电极	UU-DLC-A	2 pcs/box	B6	F
Adult suction electrodes, Universal 成人胸吸球, 通用	UU-SE-A	6 pcs/box	B8	C
Pediatric suction electrodes 儿童胸吸球	UU-SE-P	6 pcs/box	B8	D
Adult blue bulb 成人蓝色吸球	Y05-F027	6 pcs/box	A8	G
Pediatric green bulb 儿童绿色吸球	Y05-F022	6 pcs/box	A8	H





ECG Adapter 心电转接线	REF# 型号	Package 包装	Price code 价格代码	REF to Figures 图例
*Flat clip/red Adaptor/w 4mm banana 平式转换头/红色/ 4mm 香蕉头	AD-BF-10	10 pcs/box	B7	A
*Flat clip/blue adaptor/ w 3mm needle 平式转换头/蓝色/ 3mm 针头	AD-NF-10	10 pcs/box	B7	B
Alligator clip/red Adaptor/w 4mm banana 鳄鱼夹/ 红色/ 4mm 香蕉头	AD-BA-10	10 pcs/bag	B7	C
Alligator clip/blue adaptor/w 3mm needle 鳄鱼夹/ 蓝色/ 3mm 针头	AD-NA-10	10 pcs/box	B7	D

<https://endo.id>



PT. ENDO Indonesia -- Importer, distributor, and manufacturer of medical products in Indonesia

Adapter/4mm with grabber end 转接头/ 4mm 钳式	AD-BP-10	10 pcs/bag	B8	E
Adapter/4mm with snap end 转接头/ 4mm 扣式	AD-BS-10	10 pcs/box	B8	F
Adapter/4mm with banana end 转接头/ 4mm 香蕉头	AD-BW-10	10 pcs/bag	B8	I
Wide mouth Multi-function adapter/3mm with banana end 多功能心电转换钳/3mm 香蕉头	AD-NW-10	10 pcs/bag	B8	J
Adapter wire/4mm with grabber end, 0.3m 转接线/ 4mm/香蕉头 >> 钳式	AD-B-P	10 pcs/bag	B8	G
Adapter wire/3mm with grabber end, 0.3m 转接线/ 3mm/针式 >> 钳式	AD-N-P	10 pcs/box	B8	G
Adapter wire/4mm with snap end, 0.3m 转接线/ 4mm/香蕉头 > 扣式	AD-B-S	10 pcs/bag	B8	H
Adapter wire/3mm with snap end, 0.3m 转接线/ 3mm/针式 >> 扣式	AD-N-S	10 pcs/box	B8	H

Items indicated “*” are not sold in the USA.

带*的产品均不在美国区域销售

