


## Disposable Endoscopic Cutter Stapler and Cartridge

- Applied to tissue transection, resection and anastomosis in thoracic surgery, abdominal surgery and gynecologic endoscopy.



The operation is safer by resetting the button manually.



Indicate the direction of the blade walking to ensure safe operation without accident.



Adaptive joint head can swing within 90°.The joint moves flexible and can be operated by single hand conveniently.



3+1 firing model, manual firing back after 3 step of firing, safer and more stable.

### Stapler

Product Code	Model	Basic Size $\pm 10$
12018345	ECS-45-B(R)	516
12018445	ECS-45-L(R)	616
12018360	ECS-60-B(R)	528
12018460	ECS-60-L(R)	628

### Cartridge

Product Code	Model	Effective Length $\pm 2$	Color	Staple Height $\pm 0.2$	Staple Qty	Trocar
12028045	ECZ-45-M	48.5	M Grey	2.0	70	12
12028145	ECZ-45-W		W White	2.5		
12028245	ECZ-45-B		B Blue	3.5		
12028345	ECZ-45-D		D Gold	3.8		
12028445	ECZ-45-G	60	G Green	4.1	88	15
12028060	ECZ-60-W		W White	2.5		
12028460	ECZ-60-B		B Blue	3.5		
12028560	ECZ-60-D		D Gold	3.8		
12028360	ECZ-60-G		G Green	4.1		15